


INTERFACING EMCCD CAMERA USED IN ION TRAP EXPERIMENTS WITH PYTHON

YUDI WU

IMPERIAL COLLEGE LONDON

16/07/2019




- Iontraps and EMCCD Cameras
- Motivation and Aims of the Project
- Current Progress
- Future plans

Iontraps and EMCCD Cameras

- An ion trap is a device that is used to levitate small clouds of ions, or a single atomic ion, in free space, inside a vacuum chamber.
- For example a Penning trap confines the motion of the ions through the use of static electric and magnetic fields

- A Charged Coupled Device (CCD) is a silicon based semiconductor chip which captures light and converts the photons to digital data in the form of electrons
- An Electron Multiplying CCD (EMCCD) has an identical structure to conventional CCDs but is more sensitive and capable of single photo detection
- The shift register in EMCCDs are extended to include the Gain register which significantly improves low light detection

EMCCD CAMERAS



Motivation and Aims of the Project

MOTIVATION AND AIMS

- Write a program in Python to interface an Andor iXon Ultra EMCCD camera used in ion trap experiments.
- Find a method of distinguishing between a bright ion and a dark ion
- Find the setting of the program which takes the best images of the ions

Current Progress

THE PROGRAM

- Two separate programs were written in Python: one to control the EMCCD camera and take pictures of ion(s) and the other to load the images
- A Dynamic Link Library (DLL) containing various functions of the camera used by the program to control the camera
- Both programs are object oriented and have its own graphical user interface (GUI) created using QT Creator software

CAMERA CONTROL PROGRAM

cam_gui

Initialisation: OFF Camera ON Camera OFF

Temperature: -70 °C OFF -70 Set

Acquisition Mode: N/A Single Scan Frame Transfer OFF ON OFF

Read Mode: N/A Image

Exposure Time: 0 s 0.100

EMCCD Gain: 0 300

Kinetic Settings

Number of Accumulations: 1 1

Accumulation Cycle Time: 0.000 s 0.000

Number in Kinetic Series: 1 1

Kinetic Cycle Time: 0 s 0.000

Image Area

Start Column: 1 1

End Column: 512 512

Start Row: 1 1

End Row: 512 512

Trigger Mode: N/A Internal Set Trigger Mode

Set Vertical Clock Voltage Amplitude: 0 V 2 Set Voltage

EM Gain Mode: 0 0 Set EM Gain Mode

Pre-Amplifier Gain: 1 1 Set Pre-Amp Gain


Set Shutter: N/A Permanently Open Set Shutter

Vertical Shift Speed: 0.3 μ s 0.3 Set Shift Speed

Snap and Save

Snap Abort Browse Auto Save Save


IMAGE LOADING PROGRAM




INVESTIGATING CAMERA PROPERTIES

- The camera control program was first tested with noise readings and compared with commercial software to ensure the python program is working as expected
- The affect of the Exposure time and EMCCD gain on the mean noise reading per pixel were tested

MEAN SIGNAL COUNT VS. EXPOSURE TIME




(a) Python Program




(b) Commercial Software

MEAN SIGNAL COUNT VS. EMCCD GAIN




(a) Python Program




(b) Commercial Software


IMAGES OF SINGLE ION AT DIFFERENT EXPOSURE TIMES



(a) 5 Second Exposure




(b) 2 Second Exposure




(c) 1 Second Exposure


IMAGES OF SINGLE ION AT DIFFERENT EXPOSURE TIMES (CONT.)



(d) 0.7 Second Exposure



(e) 0.4 Second Exposure




(f) 0.1 Second Exposure

Future Plans

FUTURE PLANS

- Investigate the minimum exposure required for a bright ion to be able to be distinguished from a dark ion
- A bright ion can no longer be distinguished from a dark ion if the distribution of the dark ion and bright ion signal counts have a large overlap



FUTURE PLANS (CONT.)

- Comparing the differences in quality for image taken when the camera is externally triggered by the experiment and when the camera triggers the experiment and the start of exposure
- The EMCCD has a 'keep clean' cycle which clears the sensor to ensure it is charge free before the next exposure.
- Externally triggering the camera may interrupt the keep clean cycle and produce a more noisy image
- The camera gives of a 'fire signal' during exposure which can be used to trigger the start of the experiment at the end of a 'keep clean' cycle

CONCLUSIONS

- A python program with a GUI was created to acquire pictures of fluorescence of single ions and another program was created to view the images and show the vertical and horizontal projections
- A reliable method/algorithm is required to distinguish bright ions from dark ions at short exposure time
- Find the camera triggering mode which gives images of highest signal to noise ratio

Thank You for Listening!